

### Test de drogas en orina U-Cup

Prueba de detección de drogas rápida de un solo paso para la detección cualitativa de drogas y sus metabolitos en orina humana.

- CLIA waived.
- FDA approved.
- ISO13485.
- Configurable de 1 a 18 tiras reactivas.
- Opción de incluir tira de adulterantes.
- Uso rápido, los resultados aparecen en 2 a 5 minutos.
- Alta precisión.
- A prueba de filtración con mecanismo de bloqueo de tapa.
- ISO13485.
- Vida útil 24 meses.
- Almacenamiento entre 2 a 30°C

Solo para uso forense



### Interpretación de Resultados

C=Control T=Test



### Interpretación Resultados Adulterantes

TEST	ABNORMAL (LOW)	NORMAL	ABNORMAL (HIGH)
OX / PCC	Light green	White / Negative	Dark blue
S.G.	1.000 (Dark blue)	1.003, 1.005, 1.015, 1.025 (Green to Yellow)	≥1.035 (Yellow)
pH	2.0, 3.0 (Pink to Orange)	4.0, 7.0, 9.0 (Yellow to Brown)	10, 11, 12 (Dark Brown to Black)
Nitrite (mg/dl)		0, 20 (Light to Dark Orange)	50, 100 (Red to Dark Red)
GLUT	Dark Purple	Dark Purple	Dark Blue
CRE (mg/dl)	0 (Light Orange)	20, 100 (Light to Dark Purple)	